

SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006



Solution FT1 (EliGene® FFPE Tissue DNA Isolation Kit, 415050)

SECTION 1: IDENTIFICATION OF THE SUBSTANCE/PREPARATION AND THE COMPANY/UNDERTAKING

1.1 Product information

Commercial Product Name: Solution FT1 (EliGene® FFPE Tissue DNA Isolation Kit, 415050)

1.2 Relevant identified uses of the substance or mixture and uses advised against

Use of the Substance/Mixture: Laboratory chemicals

1.3 Details of the supplier of the safety data sheet

Company: ELISABETH PHARMACON, spol. s r.o.
Rokycanova 4437/5, 615 00 Brno,
Czech Republic

Telephone: +420 542 213 851

Fax: +420 542 213 827

E-mail: info@elisabeth.cz

1.4 Emergency telephone number

Centre for detoxification Prague: +420 224 91 92 93/ +420 224 91 54 02

SECTION 2: HAZARDS IDENTIFICATION

2.1 Classification of the substance or mixture

The substance or mixture is classified as dangerous according to Regulation (EC) No 1272/2008.

Acute Tox. 4: H302, H312

Skin Corr. 1B: H314

Skin Irrit. 2: H315

Aquatic Chronic Tox. 3: H412

2.2 Label elements

Hazard pictograms:



Signal Word:

Danger

Hazard Statements:

H302: Harmful if swallowed.

H312: Harmful in contact with skin.

H315: Causes skin irritation.

H319: Causes serious eye irritation.

H335: May cause respiratory irritation.

Advice on safe handling

P280 Wear protective gloves/ protective clothing/ eye protection/ face protection.

Solution FT1 (EliGene® FFPE Tissue DNA Isolation Kit, 415050)**2.3 Other hazards**

This substance/mixture contains no components considered to be either persistent, bioaccumulative and toxic (PBT), or very persistent and very bioaccumulative (vPvB) at levels of 0.1% or higher.

SECTION 3: COMPOSITION/INFORMATION ON INGREDIENTS**3.1 Hazardous ingredients**

Chemical Name	CAS- No. EC-No. Registration number	Classification (REGULATION (EC) No 1272/2008)	Concentration (%)
Guanidine thiocyanate	593-84-0 209-812-1 615-030-00-5	Acute Tox. 4: H302, H312 Skin Corr. 1B: H314 Aquatic Chronic Tox. 3: H412	>= 40 -< 50

For explanation of abbreviations see section 16.

SECTION 4: FIRST AID MEASURES**4.1 Description of first aid measures**

General advice:	Show this data sheet to the doctor in attendance.
If inhaled:	Move to fresh air. If symptoms persist, call a physician.
In case of skin contact:	Wash off immediately with soap and plenty of water while removing all contaminated clothes and shoes. If symptoms persist, call a physician.
In case of eye contact:	Remove contact lenses. Protect unharmed eye. Rinse thoroughly with plenty of water for at least 15 minutes and consult a physician.
If swallowed:	If accidentally swallowed obtain immediate medical attention. Rinse mouth with water. Never give anything by mouth to an unconscious person.

4.2 Most important symptoms and effects, both acute and delayed

Symptoms:	No information available.
-----------	---------------------------

4.3 Indication of any immediate medical attention and special treatment needed

Treatment:	No information available.
------------	---------------------------

SECTION 5: FIRE-FIGHTING MEASURES**5.1 Extinguishing media**

Suitable extinguishing media:	Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.
-------------------------------	---

5.2 Special hazards arising from the substance or mixture

Specific hazards during firefighting:	Exposure to decomposition products may be a hazard to health.
---------------------------------------	---

Solution FT1 (EliGene® FFPE Tissue DNA Isolation Kit, 415050)

Hazardous combustion products: No hazardous combustion products are known

5.3 Advice for firefighters

Special protective equipment for firefighters: Wear self-contained breathing apparatus for firefighting if necessary.

Further information: In the event of fire and/or explosion do not breathe fumes.

SECTION 6: ACCIDENTAL RELEASE MEASURES**6.1 Personal precautions, protective equipment and emergency procedures**

Personal precautions: Use personal protective equipment.
Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray.

6.2 Environmental precautions

Prevent further leakage or spillage if safe to do so.

6.3 Methods and material for containment and cleaning up

Methods for cleaning up: Keep in suitable, closed containers for disposal.

6.4 Reference to other sections

For personal protection see section 8.

SECTION 7: HANDLING AND STORAGE**7.1 Precautions for safe handling**

Advice on safe handling: Avoid contact with skin and eyes.
For personal protection see section 8.
Smoking, eating and drinking should be prohibited in the application area.

Advice on protection against fire and explosion: Normal measures for preventive fire protection.

Hygiene measures: Keep away from food and drink. Wash hands before breaks and at the end of workday. Ensure adequate ventilation, especially in confined areas. Avoid contact with the skin and the eyes. When using do not eat, drink or smoke.

7.2 Conditions for safe storage, including any incompatibilities

Requirements for storage: Keep container tightly closed in a dry and well-ventilated place.

Advice on common storage: Do not store near acids.

Other data: No decomposition if stored and applied as directed.

7.3 Specific end use(s)

Specific use(s): Laboratory chemicals

Solution FT1 (EliGene® FFPE Tissue DNA Isolation Kit, 415050)**SECTION 8: EXPOSURE CONTROLS / PERSONAL PROTECTION****8.1 Control parameters**

Contains no substances with occupational exposure limit values.

8.2 Exposure controls*Personal protective equipment*

Eye protection: Safety glasses

Hand protection Remarks: The choice of an appropriate glove does not only depend on its material but also on other quality features and is different from one producer to the other. Take note of the information given by the producer concerning permeability and break through times, and of special workplace conditions (mechanical strain, duration of contact).

Skin and body protection: Choose body protection according to the amount and concentration of the dangerous substance at the work place.
Footwear protecting against chemicals

SECTION 9: PHYSICAL AND CHEMICAL PROPERTIES**9.1 Information on basic physical and chemical properties**

Appearance:	liquid
Colour:	clear
Odour:	stench
Odour Threshold:	No data available
pH:	No data available
Melting point/range:	ca. 0 °C
Boiling point/boiling range:	ca. 100 °C
Flash point:	No data available
Evaporation rate:	No data available
Burning rate:	No data available
Upper explosion limit:	No data available
Lower explosion limit:	No data available
Vapour pressure:	No data available
Relative vapour density:	No data available
Relative density:	No data available
Density:	ca. 1 g/cm ³
Solubility(ies)	
Water solubility:	soluble
Solubility in other solvents:	No data available

Solution FT1 (EliGene® FFPE Tissue DNA Isolation Kit, 415050)

Partition coefficient: n-octanol/water: No data available

Auto-ignition temperature: No data available

Thermal decomposition: No data available

Viscosity, dynamic: No data available

Viscosity, kinematic: No data available

Explosive properties: No data available

Oxidizing properties: No data available

9.2 Other information

No data available

SECTION 10. STABILITY AND REACTIVITY**10.1 Reactivity**

No decomposition if stored and applied as directed.

10.2 Chemical stability

No decomposition if stored and applied as directed.

10.3 Possibility of hazardous reactions

Hazardous reactions: Stable under recommended storage conditions.
Hazardous decomposition products formed under fire conditions.
Keep away from oxidizing agents, and acidic or alkaline products.

10.4 Conditions to avoid

Conditions to avoid: No data available

10.5 Incompatible materials

Materials to avoid: No data available

10.6 Hazardous decomposition products

Hazardous decomposition products: No decomposition if stored and applied as directed.

SECTION 11. TOXICOLOGICAL INFORMATION**11.1 Toxicological information for the product**

Not classified based on available information.

Product:

Acute oral toxicity: Acute toxicity estimate: > 400 mg/kg
Method: Calculation method

Acute inhalation toxicity: No data available

Acute dermal toxicity: Acute toxicity estimate: > 1,000 mg/kg
Method: Calculation method

Solution FT1 (EliGene® FFPE Tissue DNA Isolation Kit, 415050)*Skin corrosion/irritation*

Not classified based on available information.

Remarks

May cause skin irritation and/or dermatitis.

Serious eye damage/eye irritation

Not classified based on available information.

*Product:**Remarks*

The product causes irritation of eyes, skin and mucous membranes.

Respiratory or skin sensitisation

Skin sensitisation: Not classified based on available information.

Respiratory sensitisation: Not classified based on available information.

Germ cell mutagenicity

Not classified based on available information.

Carcinogenicity

Not classified based on available information.

Reproductive toxicity

Not classified based on available information.

STOT - single exposure

Not classified based on available information.

STOT - single exposure

Not classified based on available information.

Aspiration toxicity

No data available

Further information

No data available

SECTION 12. ECOLOGICAL INFORMATION**12.1 Toxicity**

Toxicity to fish:

No data available

Toxicity to algae:

No data available

Toxicity to bacteria:

No data available

Components:

Guanidine thiocyanate

Toxicity to daphnia:

EC50 42 mg/l

Exposure time: 48 h

Solution FT1 (EliGene® FFPE Tissue DNA Isolation Kit, 415050)

12.2 Persistence and degradability

No data available

12.3 Bioaccumulative potential

Bioaccumulation:

No data available

12.4 Mobility in soil

No data available

12.5 Results of PBT and vPvB assessment

Assessment: This substance/mixture contains no components considered to be either persistent, bioaccumulative and toxic (PBT), or very persistent and very bioaccumulative (vPvB) at levels of 0.1% or higher.

12.6 Other adverse effects

Additional ecological information Remarks: No data available

SECTION 13. DISPOSAL CONSIDERATIONS**13.1 Waste treatment methods**

Waste from residues:

Send to a licensed waste management company.
Dispose of as hazardous waste in compliance with local and national regulations.

Contaminated packaging:

Dispose of as unused product.
Do not re-use empty containers.

SECTION 14. TRANSPORT INFORMATION**14.1 UN number**

Not regulated as a dangerous good

14.2 UN proper shipping name

Not regulated as a dangerous good

14.3 Transport hazard class(es)

Not regulated as a dangerous good

14.4 Packing group

Not regulated as a dangerous good

14.5 Environmental hazards

Not regulated as a dangerous good

14.6 Special precautions for user

For personal protection see section 8.

14.7 Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

Not applicable for product as supplied.

Solution FT1 (EliGene® FFPE Tissue DNA Isolation Kit, 415050)**SECTION 15. REGULATORY INFORMATION****15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture**

Regulation (EC) No 649/2012 of the European Parliament and the Council concerning the export and import of dangerous chemicals: Not applicable

REACH -Candidate List of Substances of Very High Concern for Authorisation (Article 59): Not applicable

Regulation (EC) No 1005/2009 on substances that deplete the ozone layer: Not applicable

Regulation (EC) No 850/2004 on persistent organic pollutants: Not applicable

Seveso II - Directive 2003/105/EC amending Council Directive 96/82/EC on the control of major-accident hazards involving dangerous substances: Not applicable

Seveso III: Directive 2012/18/EU of the European Parliament and of the Council on the control of major-accident hazards involving dangerous substances: Not applicable

15.2 Chemical Safety Assessment

No data available

SECTION 16. OTHER INFORMATION**Full text of H-Statements**

H302: Harmful if swallowed.
H312: Harmful in contact with skin.
H315: Causes skin irritation.
H319: Causes serious eye irritation.
H335: May cause respiratory irritation.

Full text of other abbreviations

Acute Tox.: Acute toxicity
Skin Corr.: Skin Corrosion
Skin Irrit.: Skin irritation
Aquatic Chronic Tox.: Aquatic Chronic Toxicity
STOT SE: Specific target organ toxicity -single exposure

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.