

**Solution FT3 (EliGene® FFPE Tissue DNA Isolation Kit, 415050)****SECTION 1: IDENTIFICATION OF THE SUBSTANCE/PREPARATION AND THE COMPANY/UNDERTAKING****1.1 Product information**

Commercial Product Name: Solution FT3 (EliGene® FFPE Tissue DNA Isolation Kit, 415050)

**1.2 Relevant identified uses of the substance or mixture and uses advised against**

Use of the Substance/Mixture: Laboratory chemicals

**1.3 Details of the supplier of the safety data sheet**

Company: ELISABETH PHARMACON, spol. s r.o.  
Rokycanova 4437/5, 615 00 Brno,  
Czech Republic

Telephone: +420 542 213 851

Fax: +420 542 213 827

E-mail: info@elisabeth.cz

**1.4 Emergency telephone number**

Centre for detoxification Prague: +420 224 91 92 93/ +420 224 91 54 02

**SECTION 2: HAZARDS IDENTIFICATION****2.1 Classification of the substance or mixture**

The substance or mixture is classified as dangerous according to Regulation (EC) No 1272/2008.

Flam. Liq. 2; H225

Eye Irrit. 2; H319

**2.2 Label elements**

Hazard pictograms:



Signal Word: Danger

Hazard Statements: H225 Highly flammable liquid and vapor.  
H319 Causes serious eye irritation.  
H336 May cause drowsiness or dizziness.

**Advice on safe handling**

P210 Keep away from heat/sparks/open flames/hot surfaces.  
No smoking.

P233 Keep container tightly closed.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses if present and easy to do – continue rinsing.

**2.3 Other hazards**

This substance/mixture contains no components considered to be either persistent, bioaccumulative and toxic (PBT), or very persistent and very bioaccumulative (vPvB) at levels of 0.1% or higher.

## SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006



### Solution FT3 (EliGene® FFPE Tissue DNA Isolation Kit, 415050)

## SECTION 3: COMPOSITION/INFORMATION ON INGREDIENTS

### 3.1 Mixtures

*Hazardous components*

Chemical Name	CAS- No. EC-No. Registration Number	Classification (REGULATION (EC) No 1272/2008)	Concentration (%)
ethanol	64-17-5 200-578-6	Flam. Liq. 2; H225 Eye Irrit. 2; H319	>= 99

For explanation of abbreviations see section 16

## SECTION 4: FIRST AID MEASURES

### 4.1 Description of first aid measures

General advice:	Move out of dangerous area Show this data sheet to the doctor in attendance.
If inhaled:	If unconscious place in recovery position and seek medical advice. If symptoms persist, call a physician.
In case of skin contact:	Wash off immediately with soap and plenty of water while removing all contaminated clothes and shoes. If symptoms persist, call a physician.
In case of eye contact:	Immediately flush eye(s) with plenty of water. Remove contact lenses. Protect unharmed eye. Rinse thoroughly with plenty of water for at least 15 minutes and consult a physician.
If swallowed:	If accidentally swallowed obtain immediate medical attention. Rinse mouth with water. Never give anything by mouth to an unconscious person.

### 4.2 Most important symptoms and effects, both acute and delayed

Symptoms:	No information available
Symptoms:	Causes serious eye irritation.

### 4.3 Indication of any immediate medical attention and special treatment needed

Treatment:	No information available.
------------	---------------------------

## SECTION 5: FIRE-FIGHTING MEASURES

### 5.1 Extinguishing media

Suitable extinguishing media:	Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide.
-------------------------------	--

**Solution FT3 (EliGene® FFPE Tissue DNA Isolation Kit, 415050)****5.2 Special hazards arising from the substance or mixture**

Specific hazards during firefighting: Do not allow run-off from fire fighting to enter drains or water courses.  
Exposure to decomposition products may be a hazard to health.

Hazardous combustion products: No hazardous combustion products are known

**5.3 Advice for firefighters**

Special protective equipment for firefighters: Wear self-contained breathing apparatus for firefighting if necessary.

Further information: In the event of fire and/or explosion do not breathe fumes.  
Use a water spray to cool fully closed containers.

**SECTION 6: ACCIDENTAL RELEASE MEASURES****6.1 Personal precautions, protective equipment and emergency procedures**

Personal precautions: Use personal protective equipment.  
Ensure adequate ventilation.  
Remove all sources of ignition.  
Evacuate personnel to safe areas.  
Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray.  
Beware of vapors accumulating to form explosive concentrations. Vapors can accumulate in low areas

**6.2 Environmental precautions**

Prevent product from entering drains.  
Prevent further leakage or spillage if safe to do so.

**6.3 Methods and material for containment and cleaning up**

Methods for cleaning up: Contain spillage, and then collect with non-combustible absorbent material, (e.g. sand, earth, diatomaceous earth, vermiculite) and place in container for disposal according to local / national regulations (see section 13).

**6.4 Reference to other sections**

For personal protection see section 8.

**SECTION 7: HANDLING AND STORAGE****7.1 Precautions for safe handling**

Advice on safe handling: Avoid formation of aerosol.  
Do not breathe vapors/dust.  
Avoid contact with skin and eyes.  
For personal protection see section 8.  
Smoking, eating and drinking should be prohibited in the application area.  
Provide sufficient air exchange and/or exhaust in work rooms.  
Open drum carefully as content may be under pressure.  
Dispose of rinse water in accordance with local and national regulations.

**Solution FT3 (EliGene® FFPE Tissue DNA Isolation Kit, 415050)**

Advice on protection against fire and explosion:

Take necessary action to avoid static electricity discharge (which might cause ignition of organic vapours). Use only explosion-proof equipment. Keep away from open flames, hot surfaces and sources of ignition

Hygiene measures:

Keep away from food and drink. Wash hands before breaks and at the end of workday. Ensure adequate ventilation, especially in confined areas. Avoid contact with the skin and the eyes. When using do not eat, drink or smoke.

**7.2 Conditions for safe storage, including any incompatibilities**

Requirements for storage areas and containers:

Keep container tightly closed in a dry and well – ventilated place.

Advice on common storage:

Do not store together with oxidizing and self -igniting products

Other data:

No decomposition if stored and applied as directed.

**7.3 Specific end use(s)**

Specific use(s):

Laboratory chemicals

**SECTION 8: EXPOSURE CONTROLS / PERSONAL PROTECTION****8.1 Control parameters**

Chemical Name	CAS-No.	Value	Control parameters	Basis
ethanol	64-17-5	TWA	1,000 ppm 1,920 mg/m3	GB EH40

**8.2 Exposure controls**

*Personal protective equipment*

Eye protection:

Safety glasses

Wear face - shield and protective suit for abnormal processing problems. Ensure that eyewash stations and safety showers are close to the workstation location.

Hand protection Remarks:

The choice of an appropriate glove does not only depend on its material but also on other quality features and is different from one producer to the other. Take note of the information given by the producer concerning permeability and break through times, and of special workplace conditions (mechanical strain, duration of contact). Protective gloves complying with EN 374.

Skin and body protection:

Choose body protection according to the amount and concentration of the dangerous substance at the work place.  
Footwear protecting against chemicals.  
Workers should wear antistatic footwear.

**SECTION 9: PHYSICAL AND CHEMICAL PROPERTIES****9.1 Information on basic physical and chemical properties**

Appearance:

liquid

**Solution FT3 (EliGene® FFPE Tissue DNA Isolation Kit, 415050)**

Colour:	colourless
Odour:	alcohol-like
Odour Threshold:	No data available
pH:	No data available
Melting point/range:	No data available
Boiling point/boiling range:	No data available
Flash point:	No data available
Evaporation rate:	No data available
Burning rate:	flammable
Upper explosion limit:	No data available
Lower explosion limit:	No data available
Vapour pressure:	No data available
Relative vapour density:	No data available
Relative density:	No data available
Density:	No data available
Solubility(ies)	
Water solubility:	soluble
Solubility in other solvents:	No data available
Partition coefficient: n-octanol/water:	No data available
Auto-ignition temperature:	No data available
Thermal decomposition:	No data available
Viscosity:	No data available
Explosive properties:	No data available
Oxidizing properties:	No data available

**9.2 Other information**

No data available

**SECTION 10. STABILITY AND REACTIVITY****10.1 Reactivity**

No decomposition if stored and applied as directed.

**10.2 Chemical stability**

No decomposition if stored and applied as directed.

**10.3 Possibility of hazardous reactions**

Hazardous reactions:	Stable under recommended storage conditions. Hazardous decomposition products formed under fire conditions. Vapors may form explosive mixture with air.
----------------------	---

## SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006



### Solution FT3 (EliGene® FFPE Tissue DNA Isolation Kit, 415050)

Keep away from oxidizing agents, and acidic or alkaline products  
none under normal use

Hazardous decomposition products: No decomposition if stored and applied as directed.

#### 10.4 Conditions to avoid

Conditions to avoid: Heat, flames and sparks.

#### 10.5 Incompatible materials

Materials to avoid: No data available

#### 10.6 Hazardous decomposition products

Hazardous decomposition products: No decomposition if stored and applied as directed.

### SECTION 11. TOXICOLOGICAL INFORMATION

#### 11.1 Toxicological information for the product

Product:

Acute oral toxicity: No data available

Acute inhalation toxicity: No data available

Acute dermal toxicity: No data available

Ingredients:

Acute oral toxicity: LD50/rat = 10,470 mg/kg

Acute inhalation toxicity: LC50/rat = 20000 ppm

Exposure time: 10 hours

*Skin corrosion/irritation*

Not classified based on available information

Product:

Remarks:

May irritate skin.

*Serious eye damage/eye irritation*

Causes serious eye irritation

*Product*

May cause irreversible eye damage.

*Ingredients*

Ethanol: eye irritation

*Respiratory or skin sensitisation*

Skin sensitisation: No data available

Respiratory sensitisation: No data available

*Germ cell mutagenicity*

No data available

*Carcinogenicity*

No data available

*Reproductive toxicity*

No data available

**Solution FT3 (EliGene® FFPE Tissue DNA Isolation Kit, 415050)**

---

*STOT - single exposure*

No data available

*Aspiration toxicity*

No data available

*Further information*

Solvents may degrease the skin.

**SECTION 12. ECOLOGICAL INFORMATION****12.1 Toxicity**

Toxicity to fish:	No data available
Toxicity to algae:	No data available
Toxicity to bacteria:	No data available

**12.2 Persistence and degradability**

No data available

**12.3 Bioaccumulative potential**

Bioaccumulation:	No data available
Partition coefficient: n	
-octanol/water:	No data available

**12.4 Mobility in soil**

No data available

**12.5 Results of PBT and vPvB assessment**

Assessment: This substance/mixture contains no components considered to be either persistent, bioaccumulative and toxic (PBT), or very persistent and very bioaccumulative (vPvB) at levels of 0.1% or higher.

**12.6 Other adverse effects**

Additional ecological information Remarks: Harmful for water organisms.

**SECTION 13. DISPOSAL CONSIDERATIONS****13.1 Waste treatment methods**

Waste from residues:	Send to a licensed waste management company. Dispose of as hazardous waste in compliance with local and national regulations.
Contaminated packaging:	Dispose of as unused product. Do not re-use empty containers.

**SECTION 14. TRANSPORT INFORMATION**

## SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006



### Solution FT3 (EliGene® FFPE Tissue DNA Isolation Kit, 415050)

#### 14.1 UN number

1170

#### 14.2 UN proper shipping name

Ethanol solution

#### 14.3 Transport hazard class(es)

3

#### 14.4 Packing group

Packing group: II  
Classification Code: F1  
Hazard Identification Number: 33  
Labels: 3  
Tunnel restriction code: (D/E)

#### 14.5 Environmental hazards

Not regulated as a dangerous good

#### 14.6 Special precautions for user

For personal protection see section 8.

#### 14.7 Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

Not applicable for product as supplied.

## SECTION 15. REGULATORY INFORMATION

#### 15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

Regulation (EC) No 649/2012 of the European Parliament and the Council concerning the export and import of dangerous chemicals: Not applicable

REACH -Candidate List of Substances of Very High Concern for Authorisation (Article 59): Not applicable

Regulation (EC) No 1005/2009 on substances that deplete the ozone layer: Not applicable

Regulation (EC) No 850/2004 on persistent organic pollutants: Not applicable

Seveso II - Directive 2003/105/EC amending Council Directive 96/82/EC on the control of major-accident hazards involving dangerous substances:

		Quantity 1	Quantity 2
7b	highly flammable	5,000 t	50,000 t

Seveso III: Directive 2012/18/EU of the European Parliament and of the Council on the control of major-accident hazards involving dangerous substances:

P5c	FLAMMABLE LIQUIDS	5,000 t	50,000 t
-----	-------------------	---------	----------

34	Petroleum products: (a) gasolines and naphthas, (b) kerosenes (including jet fuels), (c) gas oils	2,500 t	25,000 t
----	---	---------	----------



**Solution FT3 (EliGene® FFPE Tissue DNA Isolation Kit, 415050)**

---

(including diesel fuels,  
home heating oils and gas  
oil blending streams)

**15.2 Chemical Safety Assessment**

No data available

**SECTION 16. OTHER INFORMATION**

**Full text of H-Statements**

H225: Highly flammable liquid and vapour.

H319: Causes serious eye irritation.

**Full text of other abbreviations**

Eye Irrit.: Eye irritation

Flam. Liq.: Flammable liquids

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.