

# SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006



## ELiDNA™ PS Green Ultra (ED03, ED03s)

### SECTION 1: IDENTIFICATION OF THE SUBSTANCE/PREPARATION AND THE COMPANY/UNDERTAKING

#### 1.1 Product information

Commercial Product Name: ELiDNA™ PS Green Ultra (ED03, ED03s)

#### 1.2 Relevant identified uses of the substance or mixture and uses advised against

Use of the Substance/Mixture: Laboratory chemicals

#### 1.3 Details of the supplier of the safety data sheet

Company: ELISABETH PHARMACON, spol. s r.o.  
Rokycanova 4437/5, 615 00 Brno,  
Czech Republic

Telephone: +420 542 213 851

Fax: +420 542 213 827

E-mail: info@elisabeth.cz

#### 1.4 Emergency telephone number

Centre for detoxification Prague: +420 224 91 92 93/ +420 224 91 54 02

### SECTION 2: HAZARDS IDENTIFICATION

#### 2.1 Classification of the substance or mixture

Not a hazardous substance or mixture according to Regulation (EC) No. 1272/2008.

#### 2.2 Label elements

Hazard statements (H statements)

H315 Causes skin irritation

H319 Causes serious eye irritation

H335 May cause respiratory irritation

Precautionary statements (P statements)

P261 Avoid breathing mist / vapours / spray.

P280 Wear protective gloves / eye protection.

#### 2.3 Other hazards

This substance/mixture contains no components considered to be either persistent, bioaccumulative and toxic (PBT), or very persistent and very bioaccumulative (vPvB) at levels of 0.1% or higher.

### SECTION 3: COMPOSITION/INFORMATION ON INGREDIENTS

#### 3.1 Mixtures

*Hazardous components*

Chemical name	REACH Ref. No.	EC No.	Index No.	CAS No.	Content	Classification
Active ingredients of ELiDNA™ PS Green Ultra	Not available	Not available	Not available	Not available	< 1%	Not classified as hazardous
DMSO	Not available	200-664-3	Not available	67-85-5	90-100%	H315, H319, H335

**ELiDNA™ PS Green Ultra (ED03, ED03s)****SECTION 4: FIRST AID MEASURES****4.1 Description of first aid measures**

General advice:	Show this data sheet to the doctor in attendance.
If inhaled:	Move to fresh air. If symptoms persist, call a physician.
In case of skin contact:	Wash off immediately with soap and plenty of water while removing all contaminated clothes and shoes. If symptoms persist, call a physician.
In case of eye contact:	Remove contact lenses. Protect unharmed eye. Rinse thoroughly with plenty of water for at least 15 minutes and consult a physician.
If swallowed:	If accidentally swallowed obtain immediate medical attention. Rinse mouth with water. Never give anything by mouth to an unconscious person.

**4.2 Most important symptoms and effects, both acute and delayed**

Symptoms:	No information available.
-----------	---------------------------

**4.3 Indication of any immediate medical attention and special treatment needed**

Treatment:	No information available.
------------	---------------------------

**SECTION 5: FIRE-FIGHTING MEASURES****5.1 Extinguishing media**

Suitable extinguishing media:	Use extinguishing measures that are appropriate to local circumstances and the surrounding environment. The product is not combustible, adapt the measures to the environment: water spray, carbon dioxide, extinguishing foam and dry extinguishing agent.
-------------------------------	---

**5.2 Special hazards arising from the substance or mixture**

Specific hazards during firefighting:	Exposure to decomposition products may be a hazard to health.
Hazardous combustion products:	In the case of fire toxic gases – carbon monoxide, carbon dioxide, nitrogen oxides, hydrogen iodide – may be released.

**5.3 Advice for firefighters**

Special protective equipment for firefighters:	Wear self-contained breathing apparatus for firefighting if necessary.
Further information:	In the event of fire and/or explosion do not breathe fumes.

**SECTION 6: ACCIDENTAL RELEASE MEASURES****6.1 Personal precautions, protective equipment and emergency procedures**

Personal precautions:	Use personal protective equipment. Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray.
-----------------------	--

**ELiDNA™ PS Green Ultra (ED03, ED03s)****6.2 Environmental precautions**

Do not allow to enter sewers / surface waters / ground water. Remove fire residues and contaminated aqueous wastes in suitable containers and dispose of them in a controlled manner.

**6.3 Methods and material for containment and cleaning up**

Methods for cleaning up: Take up with liquid-binding agents, e.g. universal binder, treat material as described in section 13 Disposal considerations. Clean contaminated surfaces with water.

**6.4 Reference to other sections**

For personal protection see section 8.

**SECTION 7: HANDLING AND STORAGE****7.1 Precautions for safe handling**

Advice on safe handling: For personal protection see section 8.  
Smoking, eating and drinking should be prohibited in the application area.

Advice on protection against fire and explosion: Normal measures for preventive fire protection.

Hygiene measures: Keep away from food and drink. When using do not eat, drink or smoke.

**7.2 Conditions for safe storage, including any incompatibilities**

Other data: Nothing but substances of similar properties should be cumulatively stored. Cumulative storage with these substances is prohibited: pharmaceuticals, food and feed, including additives; infectious, radioactive and explosive substances; organic peroxides and very reactive oxidizing substances.

**7.3 Specific end use(s)**

Specific use(s): Laboratory chemicals

**SECTION 8: EXPOSURE CONTROLS / PERSONAL PROTECTION****8.1 Control parameters**

Contains no substances with occupational exposure limit values.

**8.2 Exposure controls***Personal protective equipment*

Eye protection: Safety glasses

Hand protection Remarks: The choice of an appropriate glove does not only depend on its material but also on other quality features and is different from one producer to the other. Take note of the information given by the producer concerning permeability and break through times, and of special workplace conditions (mechanical strain, duration of contact).

Skin and body protection: Choose body protection according to the amount and concentration of the dangerous substance at the work place.  
Footwear protecting against chemicals

**ELiDNA™ PS Green Ultra (ED03, ED03s)****SECTION 9: PHYSICAL AND CHEMICAL PROPERTIES****9.1 Information on basic physical and chemical properties**

Appearance:	liquid
Colour:	orange
Odour:	odourless
Odour Threshold:	No data available
pH:	7.0
Melting point/range:	No data available
Boiling point/boiling range:	100 °C
Flash point:	150 °C
Evaporation rate:	No data available
Burning rate:	No data available
Upper explosion limit:	No data available
Lower explosion limit:	No data available
Vapour pressure:	No data available
Relative vapour density:	No data available
Relative density:	No data available
Density:	ca 1.1 g/cm <sup>3</sup>
Solubility(ies)	
Water solubility:	soluble
Solubility in other solvents:	No data available
Partition coefficient: n-octanol/water:	No data available
Auto-ignition temperature:	not determined
Thermal decomposition:	No data available
Viscosity, dynamic:	No data available
Viscosity, kinematic:	No data available
Explosive properties:	No data available
Oxidizing properties:	No data available

**9.2 Other information**

No data available

**SECTION 10. STABILITY AND REACTIVITY****10.1 Reactivity**

No decomposition if stored and applied as directed.

**EliDNA™ PS Green Ultra (ED03, ED03s)**

---

**10.2 Chemical stability**

Chemically stable under the intended use and storage conditions.

**10.3 Possibility of hazardous reactions**

Hazardous reactions:	Stable under recommended storage conditions. Hazardous decomposition products formed under fire conditions. Stable under recommended storage conditions.
Hazardous reactions:	none under normal use
Hazardous decomposition products:	At high temperature carbon monoxide, carbon dioxide, nitrogen oxides, hydrogen iodide.

**10.4 Conditions to avoid**

Conditions to avoid:	High temperature.
----------------------	-------------------

**10.5 Incompatible materials**

Materials to avoid:	Strong acids and bases, oxidizing and reducing agents.
---------------------	--

**10.6 Hazardous decomposition products**

Hazardous decomposition products:	No decomposition if stored and applied as directed.
-----------------------------------	---

**SECTION 11. TOXICOLOGICAL INFORMATION****11.1 Toxicological information for the product**

Not classified based on available information.

Product:

Acute oral toxicity:	No data available
----------------------	-------------------

Acute inhalation toxicity:	No data available
----------------------------	-------------------

Acute dermal toxicity:	No data available
------------------------	-------------------

*Skin corrosion/irritation*

Not classified based on available information.

*Remarks*

May cause skin irritation in susceptible persons.

*Serious eye damage/eye irritation*

Not classified based on available information.

*Respiratory or skin sensitisation*

Skin sensitisation: Not classified based on available information.

Respiratory sensitisation: Not classified based on available information.

*Germ cell mutagenicity*

Not classified based on available information.

*Carcinogenicity*

Not classified based on available information.

*Reproductive toxicity*

Not classified based on available information.

*STOT - single exposure*

Not classified based on available information.

**ELiDNA™ PS Green Ultra (ED03, ED03s)**

---

*Aspiration toxicity*

Not classified based on available information.

*Further information*

No data available

**SECTION 12. ECOLOGICAL INFORMATION****12.1 Toxicity**

Toxicity to fish: No data available

Toxicity to algae: No data available

Toxicity to bacteria: No data available

**12.2 Persistence and degradability**

No data available

**12.3 Bioaccumulative potential**

Bioaccumulation: No data available

**12.4 Mobility in soil**

No data available

**12.5 Results of PBT and vPvB assessment**

Assessment: This substance/mixture contains no components considered to be either persistent, bioaccumulative and toxic (PBT), or very persistent and very bioaccumulative (vPvB) at levels of 0.1% or higher.

**12.6 Other adverse effects**

Additional ecological information Remarks: No data available

**SECTION 13. DISPOSAL CONSIDERATIONS****13.1 Waste treatment methods**

Contaminated packaging: Empty containers should be taken to an approved waste handling site for recycling or disposal.

**SECTION 14. TRANSPORT INFORMATION****14.1 UN number**

Not regulated as a dangerous good

**14.2 UN proper shipping name**

Not regulated as a dangerous good

**14.3 Transport hazard class(es)**

Not regulated as a dangerous good

**14.4 Packing group**

Not regulated as a dangerous good

**ELiDNA™ PS Green Ultra (ED03, ED03s)**

---

**14.5 Environmental hazards**

Not regulated as a dangerous good

**14.6 Special precautions for user**

For personal protection see section 8.

**14.7 Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code**

Not applicable for product as supplied.

**SECTION 15. REGULATORY INFORMATION****15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture**

REACH - Restrictions on the manufacture, placing on the market and use of certain dangerous substances, preparations and articles (Annex XVII): Not applicable

Regulation (EC) No 649/2012 of the European Parliament and the Council concerning the export and import of dangerous chemicals: Not applicable

REACH - Candidate List of Substances of Very High: Not applicable

REACH - List of substances subject to authorisation (Annex XIV): Not applicable

Regulation (EC) No 1005/2009 on substances that deplete the ozone layer: Not applicable

Regulation (EC) No 850/2004 on persistent organic pollutants: Not applicable

Seveso II - Directive 2003/105/EC amending Council Directive 96/82/EC on the control of major-accident hazards involving dangerous substances: Not applicable

**15.2 Chemical Safety Assessment**

No data available

**SECTION 16. OTHER INFORMATION**

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.