

ELISABETH PHARMACON Ltd.
Rokycanova 4437/5
615 00 Brno-Zidenice, Czech Republic

Phone: +420 542 213 851 E-mail: info@elisabeth.cz Web: www.elisabeth.cz VAT: CZ26258412



EliZyme™ Ladder 50 - 1500bp

Intended use:

For Research Use Only. Not for use in diagnostic procedures.

Storage:

Store at 25 °C for 6 months. Store at 4 °C for 12 months. Store at -20 °C for 24 months.

Product description

A combination of PCR products and restriction enzyme digested plasmids. DNA amount of each band is provided based on loading 5 μ L per well. EliZymeTM Ladder 50 - 1500bp is ready to load, no loading buffer is required.

Content

	Ref. No.	Content	Size
EliZyme™ Ladder 50 -	EZ8710	1×500 μl + 1×400 μl Loading Buffer	100 lanes
1500bp	EZ8750	5×500 μl + 1×2 ml Loading Buffer	500 lanes

	Buffer	Content	
6X Sample Loading Buffer A	6X buffer	orange G	

Specifications

PCR products and double-stranded DNA digested with appropriate restriction enzymes, are phenol extracted and equilibrated to 10 mM Tris-HCl (pH 8.0) and 10 mM EDTA.

Concentration: 112 μg/mL Recommended Load: 5 μL/well




ELISABETH PHARMACON Ltd. Rokycanova 4437/5

615 00 Brno-Zidenice, Czech Republic **Phone:** +420 542 213 851

Phone: +420 542 213 851 E-mail: info@elisabeth.cz Web: www.elisabeth.cz VAT: CZ26258412





2 % TAE agarose gel

Manufacturer:

ELISABETH PHARMACON, spol. s r. o.

Rokycanova 4437/5, Brno-Židenice 615 00

info@elisabeth.cz | www.elisabeth.cz | tel.: +420 542 213 851



Catalog number



Batch code



Use by (last day of month)



Upper limit of temperature



Manufacturer



Contains sufficient "N" tests