

			αPlastic tips compatibility																		
Brand & Model	Volume Range	Type	AP-10-N-B	AP-10-N-L-B	AP-10-N-R	AP-10-N-L-R	AP-200-Y-N-B	AP-200-N-L-B	AP-200-Y-N-R	AP-200-N-L-R	AP-200-N-L-REF	AP-1000-B-N-B	AP-1000-N-L-B	AP-1000-B-N-R	AP-1000-N-L-R	AP-1250-N-B	AP-1250-N-L-B	AP-1250-N-R	AP-1250-N-L-R	AP-1250-N-L-REF	AP-100-S-L-P-R
Eppendorf																					
Eppendorf Research Plus	0.1-2.5 µL	1 ch	✓	✓	✓	✓															
Eppendorf Research Plus	0.5-10 µL	1 ch	✓	✓	✓	✓															
Eppendorf Research Plus	2-20 µL	1 ch	✓	✓	✓	✓															
Eppendorf Research Plus	10-100 µL	1 ch					✓	✓	✓	✓	✓										✓
Eppendorf Research Plus	20-200 µL	1 ch					✓	✓	✓	✓	✓										✓
Eppendorf Research Plus	120-1200	1 ch										✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Eppendorf Research Plus	0.5-10 µL	8/12 ch	✓	✓	✓	✓															
Eppendorf Research Plus	10-100 µL	8/12 ch					✓	✓	✓	✓	✓										✓
Eppendorf Research Plus	30-300 µL	8/12 ch					✓	✓	✓	✓	✓										✓
Eppendorf Research Plus	120-1200	8/12 ch										✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Thermo																					
Thermo Finnpiquette (F1)	0.2-2 µL	1 ch	✓	✓	✓	✓															
Thermo Finnpiquette (F1)	1-10 µL	1 ch	✓	✓	✓	✓															
Thermo Finnpiquette (F1)	2-20 µL	1 ch					✓	✓	✓	✓	✓										✓
Thermo Finnpiquette (F1)	5-50 µL	1 ch					✓	✓	✓	✓	✓										✓
Thermo Finnpiquette (F1)	10-100 µL	1 ch					✓	✓	✓	✓	✓										✓
Thermo Finnpiquette (F1)	20-200 µL	1 ch					✓	✓	✓	✓	✓										✓
Thermo Finnpiquette (F1)	30-300 µL	1 ch					✓	✓	✓	✓	✓										✓
Thermo Finnpiquette (F1)	100-1000	1 ch										✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Thermo Finnpiquette (F2)	0.2-2 µL	1 ch	✓	✓	✓	✓															
Thermo Finnpiquette (F2)	1-10 µL	1 ch	✓	✓	✓	✓															
Thermo Finnpiquette (F2)	2-20 µL	1 ch					✓	✓	✓	✓	✓										✓
Thermo Finnpiquette (F2)	5-50 µL	1 ch					✓	✓	✓	✓	✓										✓
Thermo Finnpiquette (F2)	10-100 µL	1 ch					✓	✓	✓	✓	✓										✓
Thermo Finnpiquette (F2)	20-200 µL	1 ch					✓	✓	✓	✓	✓										✓
Thermo Finnpiquette (F2)	30-300 µL	1 ch					✓	✓	✓	✓	✓										✓
Thermo Finnpiquette (F2)	100-1000	1 ch										✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Thermo Finnpiquette (F3)	0.2-2 µL	1 ch	✓	✓	✓	✓															
Thermo Finnpiquette (F3)	1-10 µL	1 ch	✓	✓	✓	✓															
Thermo Finnpiquette (F3)	2-20 µL	1 ch					✓	✓	✓	✓	✓										✓
Thermo Finnpiquette (F3)	5-50 µL	1 ch					✓	✓	✓	✓	✓										✓
Thermo Finnpiquette (F3)	10-100 µL	1 ch					✓	✓	✓	✓	✓										✓
Thermo Finnpiquette (F3)	20-200 µL	1 ch					✓	✓	✓	✓	✓										✓
Thermo Finnpiquette (F3)	30-300 µL	1 ch					✓	✓	✓	✓	✓										✓
Thermo Finnpiquette (F3)	100-1000	1 ch										✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Thermo Finnpiquette (F3)	0.5-10 µL	8/12 ch	✓	✓	✓	✓															
Thermo Finnpiquette (F3)	20-200 µL	8/12 ch					✓	✓	✓	✓	✓										✓
Thermo Finnpiquette (F3)	20-200 µL	8/12 ch					✓	✓	✓	✓	✓										✓
Sartorius																					
Proline Plus	0.5-10 µL	1 ch	✓	✓	✓	✓															
Proline Plus	20-200 µL	1 ch					✓	✓	✓	✓	✓										✓
Proline Plus	100-1000	1 ch										✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Proline	50-300 µL	8/12 ch					✓	✓	✓	✓	✓										✓
Gilson																					
Gilson Pipetman	0.2-2 µL	1 ch																			
Gilson Pipetman	1-10 µL	1 ch	✓	✓	✓	✓															
Gilson Pipetman	2-20 µL	1 ch																			
Gilson Pipetman	10-100 µL	1 ch					✓	✓	✓	✓	✓										✓
Gilson Pipetman	20-200 µL	1 ch					✓	✓	✓	✓	✓										✓
Gilson Pipetman	100-1000	1 ch										✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Brand																					
Transferpette S	0.1-2.5 µL	1 ch	✓	✓	✓	✓															
Transferpette S	0.5-10 µL	1 ch	✓	✓	✓	✓															
Transferpette S	2-20 µL	1 ch					✓	✓	✓	✓	✓										✓
Transferpette S	5-50 µL	1 ch					✓	✓	✓	✓	✓										✓
Transferpette S	10-100 µL	1 ch					✓	✓	✓	✓	✓										✓
Transferpette S	20-200 µL	1 ch					✓	✓	✓	✓	✓										✓
Transferpette S	30-300 µL	1 ch					✓	✓	✓	✓	✓										✓
Transferpette S	100-1000	1 ch										✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Biohit																					
Biohit Proline	0.1-2.5 µL	1 ch	✓	✓	✓	✓															
Biohit Proline	0.5-10 µL	1 ch	✓	✓	✓	✓															
Biohit Proline	2-20 µL	1 ch					✓	✓	✓	✓	✓										✓
Biohit Proline	5-50 µL	1 ch					✓	✓	✓	✓	✓										✓
Biohit Proline	10-100 µL	1 ch					✓	✓	✓	✓	✓										✓
Biohit Proline	20-200 µL	1 ch					✓	✓	✓	✓	✓										✓
Biohit Proline	50-300 µL	1 ch					✓	✓	✓	✓	✓										✓
Biohit Proline	100-1000	1 ch										✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Biohit Proline (Plus)	0.5-10 µL	1 ch	✓	✓	✓	✓															
Biohit Proline (Plus)	2-20 µL	1 ch					✓	✓	✓	✓	✓										✓
Biohit Proline (Plus)	5-50 µL	1 ch					✓	✓	✓	✓	✓										✓
Biohit Proline (Plus)	10-100 µL	1 ch					✓	✓	✓	✓	✓										✓
Biohit Proline (Plus)	20-200 µL	1 ch					✓	✓	✓	✓	✓										✓
Biohit Proline (Plus)	30-300 µL	1 ch					✓	✓	✓	✓	✓										✓
Biohit Proline (Plus)	100-1000	1 ch										✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
D-LAB																					
Micro Pette	0.1-2.5 µL	1 ch	✓	✓	✓	✓															
Micro Pette	0.5-10 µL	1 ch	✓	✓	✓	✓															
Micro Pette	2-20 µL	1 ch					✓	✓	✓	✓	✓										✓
Micro Pette	5-50 µL	1 ch					✓	✓	✓	✓	✓										✓
Micro Pette	10-100 µL	1 ch					✓	✓	✓	✓	✓										✓
Micro Pette	20-200 µL	1 ch					✓	✓	✓	✓	✓										✓
Micro Pette	50-200 µL	1 ch					✓	✓	✓	✓	✓										✓
Micro Pette	100-1000	1 ch										✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Micro Pette	200-1000	1 ch										✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
MicroPette	0.5-10 µL	8/12 ch	✓	✓	✓	✓															
MicroPette	5-50 µL	8/12 ch					✓	✓	✓	✓											