

			αPlastic tips compatibility																				
Brand & Model	Volume Range	Type	AP-10-N-B	AP-10-N-L-B	AP-10-N-R	AP-10-N-L-R	AP-200-Y-N-B	AP-200-N-L-B	AP-200-Y-N-R	AP-200-N-L-R	AP-200-N-L-REF	AP-1000-B-N-B	AP-1000-N-L-B	AP-1000-B-N-R	AP-1000-N-L-R	AP-1250-N-B	AP-1250-N-L-B	AP-1250-N-R	AP-1250-N-L-R	AP-1 250-N-L-REF	AP-100-S-L-F-R		
ZEPHYRUS																							
Zephyrus mPipet™	0.1–2.5 µL	1 ch	✓	✓	✓	✓																	
Zephyrus mPipet™	0.5–10 µL	1 ch	✓	✓	✓	✓																	
Zephyrus mPipet™	2–20 µL	1 ch																					
Zephyrus mPipet™	5–50 µL	1 ch																					
Zephyrus mPipet™	10–100 µL	1 ch					✓	✓	✓	✓											✓		
Zephyrus mPipet™	20–200 µL	1 ch					✓	✓	✓	✓	✓												
Zephyrus mPipet™	100–1000	1 ch										✓	✓	✓	✓	✓	✓	✓	✓	✓			
Zephyrus mPipet™	0.5–5 mL	1 ch																					
Zephyrus mPipet™	1–10 mL	1 ch																					
Eppendorf																							
Eppendorf Research Plus	0.1–2.5 µL	1 ch	✓	✓	✓	✓																	
Eppendorf Research Plus	0.5–10 µL	1 ch	✓	✓	✓	✓																	
Eppendorf Research Plus	2–20 µL	1 ch	✓	✓	✓	✓																	
Eppendorf Research Plus	10–100 µL	1 ch					✓	✓	✓	✓	✓										✓		
Eppendorf Research Plus	20–200 µL	1 ch					✓	✓	✓	✓	✓												
Eppendorf Research Plus	120–1200	1 ch										✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		
Eppendorf Research Plus	0.5–10 µL	8/12 ch	✓	✓	✓	✓																	
Eppendorf Research Plus	10–100 µL	8/12 ch					✓	✓	✓	✓	✓										✓		
Eppendorf Research Plus	30–300 µL	8/12 ch					✓	✓	✓	✓	✓										✓		
Eppendorf Research Plus	120–1200	8/12 ch										✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		
Thermo																							
Thermo Finnpiipette (F1)	0.2–2 µL	1 ch	✓	✓	✓	✓																	
Thermo Finnpiipette (F1)	1–10 µL	1 ch	✓	✓	✓	✓																	
Thermo Finnpiipette (F1)	2–20 µL	1 ch					✓	✓	✓	✓	✓										✓		
Thermo Finnpiipette (F1)	5–50 µL	1 ch					✓	✓	✓	✓	✓										✓		
Thermo Finnpiipette (F1)	10–100 µL	1 ch					✓	✓	✓	✓	✓										✓		
Thermo Finnpiipette (F1)	20–200 µL	1 ch					✓	✓	✓	✓	✓										✓		
Thermo Finnpiipette (F1)	30–300 µL	1 ch					✓	✓	✓	✓	✓										✓		
Thermo Finnpiipette (F1)	100–1000	1 ch										✓	✓		✓	✓	✓	✓	✓	✓	✓		
Thermo Finnpiipette (F2)	0.2–2 µL	1 ch	✓	✓	✓	✓																	
Thermo Finnpiipette (F2)	1–10 µL	1 ch	✓	✓	✓	✓																	
Thermo Finnpiipette (F2)	2–20 µL	1 ch					✓	✓	✓	✓	✓										✓		
Thermo Finnpiipette (F2)	5–50 µL	1 ch					✓	✓	✓	✓	✓												

[illegible]