

αPlastic Pipette Tips COMPATIBILITY CHART

Pipettor Brand	Pipettor Make	Channels	Volume Range	0.5-10 µL NT-10XN-S-L-F-R	0.5-10 µL C-10-S-L-F-R	1-200 µL NT-200-S-L-F-R	1-200 µL C-200-S-L-F-R	100-1 000 µL NT-1000-S-L-F-R	100-1 000 µL C-1250-S-L-F-R
Axygen®	Axyjet®	Single-channel	0.1-2 µL		✓				
			0.5-10 µL		✓				
			0.2-2 µL						
			5-50 µL						
			10-100 µL						
			20-200 µL				✓		
			50-250 µL						
			100-1,000 µL						
		Multi-channel	1-10 µL		✓				
			5-50 µL						
			20-200 µL				✓		
			50-300 µL						
	Axyjet-16	16-channel	0.2-2 µL						
			0.5-10 µL						
			5-50 µL						
Sartorius / Biohit	Proline®	Single-channel	0.1-2.5 µL	✓	✓				
			0.2-10 µL	✓	✓				
			0.5-10 µL	✓	✓				
			5-50 µL				✓		
			5-100 µL						
			10-250 µL				✓		
			10-500 µL						
			20-200 µL			✓	✓		
			50-1,000 µL			✓		✓	
			100-1000 µL					✓	
		Multi-channel	0.2-10 µL		✓				
			5-100 µL						
			25-250 µL				✓		
			50-250 µL				✓		
			5-100 µL						

Pipettor Brand	Pipettor Make	Channels	Volume Range	0.5-10 µL NT-10XN-S-L-F-R	0.5-10 µL C-10-S-L-F-R	1-200 µL NT-200-S-L-F-R	1-200 µL C-200-S-L-F-R	100-1 000 µL NT-1000-S-L-F-R	100-1 000 µL C-1250-S-L-F-R
Sartorius / Biohit	eLINE®	Single-channel	0.1-10 µL	✓	✓				
			0.2-10 µL	✓	✓				
			5-120 µL	✓		✓			
			10-300 µL						
			50-1,000 µL					✓	
		Multi-channel	0.2-10 µL		✓				
			5-120 µL						
			10-300 µL						
	mLINE®	Single-channel	0.1-3 µL		✓				
			0.5-10 µL		✓				
			2-20 µL						
			10-100 µL						
			20-200 µL			✓	✓		
			100-1,000 µL					✓	
		Multi-channel	0.5-10 µL		✓				
			5-100 µL						
			30-300 µL				✓		
BrandTech®	Transferpette®	Single-channel	0.1-1 µL	✓	✓				
			0.5-10 µL		✓				
			2-20 µL 'Ultra Micro'						
			2-20 µL						
			5-50 µL						
			10-100 µL						
			20-200 µL			✓	✓		
			25-250 µL				✓		
			100-1,000 µL					✓	
		Multi-channel	0.5-10 µL		✓				
			2-20 µL		✓				
			2.5-25 µL		✓				
			5-50 µL						
			10-100 µL						
			20-200 µL						
			30-300 µL				✓		

Pipettor Brand	Pipettor Make	Channels	Volume Range	0.5-10 µL NT-10XN-S-L-F-R	0.5-10 µL C-10-S-L-F-R	1-200 µL NT-200-S-L-F-R	1-200 µL C-200-S-L-F-R	100-1 000 µL NT-1000-S-L-F-R	100-1 000 µL C-1250-S-L-F-R
Capp Aero	Capp Aero	Single-channel	0.2-2 µL		✓				
			0.5-10 µL		✓				
			5-50 µL						
			10-100 µL						
			25-200 µL			✓			
			50-300 µL				✓		
			100-1,000 µL				✓	✓	
	Capp Aero 384 format 16/48/64	Multi-channel	0.2-2 µL						
			0.5-10 µL						
			5-50 µL						
Eppendorf®	Electronic	Single/Multi-channel	0.5-10 µL		✓				
			5-100 µL						
			2-200 µL				✓		
			50-1,000 µL						
	Reference	Single/Multi-channel	0.1-2.5 µL		✓				
			0.1-10 µL	✓	✓				
			0.5-10 µL		✓				
			2-20 µL “UltraMicro”						
			2-20 µL						
			10-100 µL						
			20-200 µL				✓		
			50-200 µL			✓	✓		
			100-1,000 µL					✓	
	Research	Single-channel	0.1-2.5 µL	✓	✓				
			0.5-10 µL	✓	✓				
			2-20 µL			✓			
			10-100 µL			✓			
			20-200 µL			✓	✓		
			100-1,000 µL					✓	
		Multi-channel	0.5-10 µL	✓	✓				
			5-100 µL			✓			
			10-100 µL			✓			
			50-100 µL				✓		

Pipettor Brand	Pipettor Make	Channels	Volume Range	0.5-10 µL NT-10XN-S-L-F-R	0.5-10 µL C-10-S-L-F-R	1-200 µL NT-200-S-L-F-R	1-200 µL C-200-S-L-F-R	100-1 000 µL NT-1000-S-L-F-R	100-1 000 µL C-1250-S-L-F-R
Eppendorf®	Research	Multi-channel	20-300 µL			✓	✓		
			30-300 µL			✓	✓		
			50-1,200 µL					✓	
	Research Pro	s/ch	20-300 µL				✓		
			50-1,000 µL					✓	
		m/ch	0.5-10 µL		✓				
			20-300 µL				✓		
	Research Plus	s/ch	100-1,000 µL					✓	
Thermo Scientific	Finnpipette™	Single-channel	0.2-2 µL		✓				
			0.3-3 µL		✓				
			0.2-10 µL		✓				
			0.5-5 µL		✓				
			0.5-10 µL		✓				
			1-10 µL		✓				
			5- 40 µL						
			5- 50 µL						
			10- 100 µL						
			40-200 µL				✓		
			30-300 µL				✓		
			100-1,000 µL					✓	
			200-1,000 µL						
	Multi-channel		1-10 µL		✓				
			5- 50 µL						
			50-300 µL				✓		
	Finnpipette Digitl	s/ch	40-200 µL				✓		
	Finnpipette electronic	Multi-channel	1-10 µL		✓				
			30-300 µL				✓		
			100-1,000 µL					✓	
			100-1200 µL						
Gilson	Pipetman® Classic	Single-channel	P2, 0.1-2 µL		✓				
			P2, 0.2-2 µL		✓				
			P10, 0.2-10 µL	✓	✓				

Pipettor Brand	Pipettor Make	Channels	Volume Range	0.5-10 µL NT-10XN-S-L-F-R	0.5-10 µL C-10-S-L-F-R	1-200 µL NT-200-S-L-F-R	1-200 µL C-200-S-L-F-R	100-1 000 µL NT-1000-S-L-F-R	100-1 000 µL C-1250-S-L-F-R
Gilson	Pipetman® Classic	Single-channel	P20, 2-20 µL						
			P100, 10-100 µL						
			P200, 20-200 µL			✓	✓		
			P1,000, 100-1,000 µL					✓	
	Pipetman Neo	Multi-channel	P2N, 0.2-2 µL		✓				
			P10N, 1-10 µL		✓				
			P20N, 2-20 µL						
			P100N, 10-100 µL						
			P100N, 20-100 µL						
			P200N, 20-200 µL				✓		
	Pipetman Neo	m/ch	P8x20N, 2-20 µL				✓		
	Pipetman Ultra	s/ch	U200, 20-200 µL				✓		
		m/ch	P200-M8, 20-200 µL		✓		✓		
Labnet	Labpette™ Discovery or	Single-channel	0.1-2 µL	✓	✓				
			0.5-10 µL		✓				
			2-20 µL						
			10-100 µL						
			20-200 µL			✓	✓		
			100-1,000 µL					✓	
		Multi-channel	1-10 µL		✓				
			5-50 µL						
			20-200 µL				✓		
			50-300 µL				✓		
	BioPette E	Single-channel	0.5-10 µL		✓				
			2-20 µL						
			10-200 µL				✓		
			100-1,000 µL						
		Multi-channel	0.5-10 µL		✓				
			2-20 µL						
			10-200 µL				✓		
			100-1200 µL						

Pipettor Brand	Pipettor Make	Channels	Volume Range	0.5-10 µL NT-10XN-S-L-F-R	0.5-10 µL C-10-S-L-F-R	1-200 µL NT-200-S-L-F-R	1-200 µL C-200-S-L-F-R	100-1 000 µL NT-1000-S-L-F-R	100-1 000 µL C-1250-S-L-F-R
Matrix	Impact / Impact2	Single/Multi-chnl	2-125 µL						
			5-250 µL				✓		
			15-1250 µL						
Nichiryo/Oxford	Benchmate	Single-channel	0.5-10 µL	✓	✓				
			2-20 µL						
			10-50 µL						
			10-100 µL						
			40-200 µL			✓			
			100-1,000 µL					✓	
		m/ch	50-300 µL						
Rainin	EDP3	Single-channel	E-MIC-10, 0.5-10 µL		✓				
			E2-25, 2.5-25 µL						
			E2-100, 10-100 µL						
			E2-250, 25-250 µL				✓		
			E2-1,000, 100-1,000 µL					✓	
		m/ch	EP-M8-10, 0.5-10 µL		✓				
			EP-M8-250, 10-250 µL				✓		
	Latch Mode	Single-channel	R2, 0.1-2 µL		✓				
			R10, 0.2-10 µL		✓				
			R20, 2-20 µL			✓			
			R100, 10-100 µL						
			R200, 20-200 µL				✓		
			R1,000, 100-1,000 µL					✓	
Socorex	Calibra/Accura	Single-channel	0.1-2 µL		✓				
			0.2-2 µL		✓				
			0.5-10 µL		✓				
			1-10 µL		✓				
			2-20 µL						
			10-100 µL						
			20-200 µL			✓			
			100-1,000 µL				✓	✓	
		Multi-channel	1-10 µL		✓				
			5-50 µL						

Pipettor Brand	Pipettor Make	Channels	Volume Range	0.5-10 µL NT-10XN-S-L-F-R	0.5-10 µL C-10-S-L-F-R	1-200 µL NT-200-S-L-F-R	1-200 µL C-200-S-L-F-R	100-1 000 µL NT-1000-S-L-F-R	100-1 000 µL C-1250-S-L-F-R
Socorex	Calibra/Accura	Multi-channel	10-100 µL						
			20-200 µL				✓		
			50-200 µL				✓		
VWR	VWR and VWR Signature	Single-channel	0.1-2 µL		✓				
			0.5-10 µL		✓				
			2-20 µL						
			5-50 µL						
			10-100 µL						
			20-200 µL				✓		
			50-250 µL				✓		
			100-1,000 µL						
		Multi-channel	1-10 µL		✓				
			5-50 µL						
			20-200 µL				✓		
			50-300 µL				✓		

αPlasticTM
Labware ELISABETH PHARMACON



ELISABETH PHARMACON, spol. s r.o.
Rokycanova 4437/5
615 00 Brno-Židenice, Česká republika
IČ: 26258412, DIČ: CZ26258412
Tel.: +420 542 213 851
E-mail: info@elisabeth.cz
www.elisabeth.cz

Kontakt:
obchod@elisabeth.cz
tel. +420 542 213 851